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Specifications

Specifications

Electrical requirements	110V 60Hz single phase, 2300 watts, 21A or 220V 50/60Hz single phase, 2300 watts, 10.5A Rotation motor: 300 watts
Glass material	High borosilicate glass 3.3
Evaporating flask	Capacity: 1.32 gallon (5 liters) Maximum load: Not to exceed 60% of the flask's capacity Sealing flange: 2" diameter Rotation speed: 10 to 140 rpm Evaporating speed: 0.5 gallon/hour (water), 1 gallon/hour (alcohol)
Receiving flask	0.79 gallon (3 liter)
Condenser	Glassware: Vertical cold trap (dry ice) condenser with lid Vacuum connection: 3/8" hose barb x 1
Working environment	5 to 95F
Water bath	Temperature range: ambient to 99C (210F) Temperature controller: digital Dimensions: 11.4" diameter x 6.7" depth x 9" height Material: stainless steel Media: distilled water only Max. water level: not to exceed 1" from the top rim of the bath Number of heaters: one

	Drain port: one Exterior liner: rubber Lifting: motoriezed Maximum lifting distance: 4 inches
Vacuum	Vacuum level: < 3 torrs Sealing: PTFE + Viton gaskets Connection: 3/8" hose barb x 1
Safety	Dry-run, over-temp shut-off, built-in circuit breakers
Shipping weight	105 Lb (evaporator), 57 Lb (glassware)
Unit dimensions	30 x 15.7 x 42" (WxDxH)
Shipping dimensions	29 x 17 x 35" (evaporator), 22 x 22 x 22" (glassware) (WxDxH)
Compliance	UL 61010-1: 2012 Ed.3 + R: 29 Apr 2016 CSA C22.2# 61010-1-12: 2012 Ed.3+ U1; U2 (R2017) CE

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A rotary evaporator is a device used in chemical laboratories for the efficient and gentle removal of solvents from samples by evaporation. The process of rotary evaporation is most often used to separate solvents with low boiling points, such as n-hexane or ethyl acetate, from compounds which are solid at room temperature and pressure. However, careful application also allows for the removal of a solvent from a sample containing a liquid compound, if there is minimal co-evaporation (azeotropic behavior) and a sufficient difference in boiling points at the chosen temperature and reduced pressure.

Example: The following table illustrates a variety of substances, ordered by increasing vapor pressure (in absolute units).

Rotary Evaporator Frequently Asked Questions

- **Q: What type of plug is needed for my rotary evaporator?**
A: The type of plug will depend on the model of rotary evaporator and whether the unit is a 110V or 220V model. Please consult the table below for further details:

Rotary Evaporator Model	
Voltage (V)	
Plug Type	
2L	
110/220	
Bath: 6-15P (220V) Motor Housing: 5-15P (110V)	
5L	
110	
L5-30P	
220	
6-20P	
10L	
110	

not available	
220	
L6-30P	
20L	
110	
not available	
220	
L6-30P	
50L	
110	
not available	
220	
L6-30P	

• **Q: What is the temperature range of the rotary evaporator**

A: The temperature range of the rotary evaporator is 20-99C.

• **Q: What is the maximum pressure that the glass condenser's coil is rated to? What is the maximum pressure that the condenser itself is rated to?**

A: The glass condenser's coil is rated to 14 psi and the condenser itself is only rated for atmospheric pressure and a vacuum of 3 torr.

• **Q: What is the dual receiving flask kit used for on the 20 and 50L rotary evaporators?**

A: The dual receiving flask kit is used to remove reclaimed solvent from the evaporator without interrupting its operation.

• **Q: What heating medium can be used in the water bath?**

A: Only distilled water can be used as a heating medium in the water bath of an Across International rotary evaporator.

• **Q: Can the vacuum level be set on the rotary evaporator?**

A: The vacuum level can be set only when a vacuum controller is purchased with the rotary evaporator. Please speak to sales or check the rotary evaporator's customization section on the website for recommendations.

• **Q: What accessories come with the purchase of a rotary evaporator?**

A: The accessories that come with each model will be different depending on the unit purchased. An accessories list can be found on the unit's website listing and can be found in the standard package tab.

• **Q: What is the evaporation speed of my rotary evaporator?**

A: The evaporation speed of the rotary evaporator will depend on what is being evaporated, and the parameters chosen. Please refer to the table for the estimated evaporation rate for ethanol and water.

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\$3,000

Certificaciones

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Across International SolventVap Rotary Evaporators, UL/CSA compliant

Ai 5.3 Gallon / 20L SolventVap RotoVap SetUp